



System 50 - Aluminium (AL)

- Round hood recommended for light fumes, gases and small open vessels
- Diameter of hood: $\varnothing 280$ mm
- Increased stability when moving the hood due to reinforced rim of the hood
- Increased efficiency at an angled position
- Hood made of transparent PETG, resistant to solvents
- Connection tube made of anodized aluminium - surface 10 μ
- The transparent hood assists in keeping a good view of the work place
- Colours of flange available: white (-5), red (-4) or black (-050)



Capture efficiency

Measurements for this product are comparable with article no. [1-6328](#).

Drawing

