



## System 75 - Aluminium (AL)

- Round hood recommended for light fumes, gases and over small open vessels
- Diameter of hood: Ø280 mm
- Increased stability when moving the hood due to reinforced rim of the hood
- Increased efficiency when placed at an angled position
- Hood made of transparent PETG, resistant to solvents
- Connection tube made of anodized aluminium - surface 10 µ
- The transparent hood assists in keeping a good view of the work place
- Colours of flange available: white (-5) or red (-4)



## Capture efficiency

Measurements for this product are comparable with article no. [1-6328](#).

## Drawing

