



## System 25

All components RoHS-compatible according to directive 2002/95/EC (RoHS).

### Design:

- Work cabinet for placing on a table/plane surface
- Dimensions of the base: 1000 x 6000 mm
- Height of side panels: 200, 500 or 750mm (total heights: 410, 710, 960 mm)
- Vacuum moulded top made of transparent PETG with a pane of tempered glass
- Pane corners made of PP in the colours white (-5) or red (-4)
- Sides made of transparent PETG
- Corners made of anodized aluminium
- Bottom edge provided with a non-skid rubber edge
- Recommended airflow: 220, 550 or 850 m<sup>3</sup>/h dependent on the height of the side panels
- Optional flange for extraction available
- NB! NB! Can be delivered with doors



## Working area

