



System 25

All components RoHS-compatible according to directive 2002/95/EC (RoHS).

Design:

- Work cabinet for placing on a table/plane surface
- Dimensions of the base: 420 x 320 mm
- Height of side panels: 400 total height: 515 mm)
- Vacuum moulded top made of transparent PETG inclusive of 2 stub pipes
- Sides made of steel with a white polyester powder-coating on the surface
- Flanges and joint with damper in the colours white (-5) or red (-4)
- Flanges, stub pipes and joint with damper $\varnothing 75$ mm
- Recommended airflow: 204 m³/h



Working area

