





Flat hood kit

Design:

- Flat hood with bracket recommended when extracting light and heavy gases and fumes
- Dimension of the hood: see table below
- Hood made of transparent PETG or black polypropylene
- Connection tube made of anodized aluminium
- Bracket made of steel with a black polyester powder-coating
- For placement on a table or hanging on a wall



Article numbers

| Alsident® System Flat hood kit | | | | |
|---|---------------------------------|---------------------------------|---------------------------------|----------------------------------|
| Transparent | System 50 (330x240mm) | System 63 (330x240mm) | System 75 (330x240mm) | System 100 (420x280mm) |
|  | 1-50-33-5 | 1-63-33-5 | 1-75-33-5 | 1-100-42-5 |
| black | System 50 | | System 75 | System 100 |
|  | 1-50-33 | | 1-75-33 | 1-100-42 |